

Specifications MicroSwing® 6.0

Manufacturer	Haynl-Elektronik GmbH, Germany
Distributor	HAIDER BIOSWING GmbH, Germany
Initial marketing	August 2012 (MicroSwing 5.0: 1999)
Type designation MicroSwing®	6.0
Serial number	ongoing
Medical product in Europe	yes
Medical product classification	1
Identification	CE
Disinfected	Wipe disinfection according to VAH Germany
Sterilizable	no
Calibration required	no
Mandatory verification	no
Type designation measuring box	HWB11
External dimensions measuring box (h x w x d)	27 x 65 x 85 mm
Mass measuring box	82 g
Type designation sensor	HCS11
sensor parameters	acceleration
recording dimensions	3 (x-, y- und z-axis)
sensor bond	cable
External dimensions sensor (h x w x d)	14 x 21 x 29 mm
Mass sensor	60 g (incl. 3,3 m cable and plug)
Accessories	Software-CD, sensor bracket, cable Protection
Medical Application	sensorimotor diagnostic assessment

BIOSWING Therapiesysteme