

Specifications BIOSWING® Posturomed® compact

Manufacturer / Distributor	HAIDER BIOSWING GmbH, Germany
Initial marketing	1997
Type designation	compact
Serial number	ongoing
Medical product in Europe	yes
Medical product classification	1
Identification	CE
Disinfected	Wipe disinfection according to VAH Germany
Sterilizable	no
Calibration required	no
Mandatory verification	no
Security barrier	one-sided folding, flexible stored
Transport wheels	2 piece / Ø 100 mm
External dimensions (unfolded, h x w x d)	994 x 699 x 614 mm
External dimensions (folded, h x w x d)	285 x 800 x 614 mm
Surface treatment (w x d)	400 x 400 mm
Maximum deflection of the surface x-axis	25 mm (locked) / 50 mm (unlocked) from the middle
Maximum deflection of the surface y-axis	40 mm (locked) / 80 mm (unlocked) from the middle
Force curve of the deflection of the surface	from the middle of exponentially growing
Oscillation frequency spectrum locked	3.8 to 6.0 Hz
Oscillation frequency spectrum unlocked	2.0 to 4.2 Hz
Total weight	26.8 kg
Maximum load capacity of the surface	180 kg
Accessories	various additional modules
Medical Application	sensorimotor therapy

BIOSWING Therapiesysteme