

## Specifications BIOSWING® Propriomed 1, Propriomed 2, Propriomed 100

Manufacturer / Distributor	HAIDER BIOSWING GmbH, Germany
Initial marketing	1999
Type designation	1, 2 and 100
Serial number	ongoing
Medical product in Europe	yes
Medical product classification	1
Identification	CE
Certification	TÜV Thüringen (Germany) international certification safety in DIN 32935, PAK (German industry standard)
Disinfected	Wipe disinfection according to VAH Germany
Sterilizable	no
Calibration required	no
Mandatory verification	no
Lengths of the vibrating wires	Propriomed 1: 190 cm, Propriomed 2: 165 cm Propriomed 100: 100 cm
Maximum amplitudes	Propriomed 1: 60 cm, Propriomed 2: 60 cm Propriomed 100: 40 cm
Spectrum of vibration frequencies	Propriomed 1: infinitely variable 3.0 to 4.0 Hz Propriomed 2: infinitely variable 4.0 to 5.2 Hz Propriomed 100: infinitely variable 4.8 to 6.0 Hz
Total weight	Propriomed 1: 975 g, Propriomed 2: 920 g Propriomed 100: 550 g
Accessories	Devices for hanging storage
Medical Application	sensorimotor therapy

BIOSWING Therapiesysteme