

Specifications BIOSWING® Torsiomed®

Manufacturer / Distributor	HAIDER BIOSWING GmbH, Germany
Initial marketing	2004
Type designation	Torsiomed
Serial number	ongoing
Medical product in Europe	yes
Medical product classification	1
Identification	CE
Disinfected	Wipe disinfection according to VAH Germany
Sterilizable	no
Calibration required	no
Mandatory verification	no
Security barrier	one-sided frontal
Transport wheels	2 piece / Ø 75 mm
External dimensions (h x w x d)	1127 x 530 x 770 mm
Surface treatment Ø	320 mm
Maximum deflection of the surface	18 mm from the middle / 36 mm total
Maximum torsion of the surface	25° from the middle / 50° total
Force curve of the deflection of the surface	from the middle of exponentially growing
Oscillation frequency spectrum	1.5 to 6.0 Hz
Total weight	23.4 kg
Maximum load capacity of the surface	150 kg
Accessories	additional modules
Medical Application	sensorimotor therapy